

RP11-191L9.1-FLJ46257 protein Gene

Reference: PTXBC132762

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FLJ46257 protein
GENE ID	400932
ORIGIN SPECIES	Human
ACCESSIONS	BC132762
VECTOR	pDONR223
NCBI SYMBOL	RP11-191L9.1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccttgcaggtttcaccaagcaagagcggagaagacacagcctggacaaaggaggtggacgggacagcagcctgcagcagcctcacaggacacagctaccagacagagagcaccagctgacactcattcctcagcaccaggaaatgccagcagaagagagtctt
caagatggccagacacacttctgcagcatggttctggagcctgcacaggccaggtttctctacaggtgcaattgagaactcactgaggcaaaaggcccttgtcaaccacagaagtgacaagtccaagtgcagaccaggttcagggtcagttctctcagaactcaaatccc
atctctgcaggtctcagccaatcctgctcccgcctgttttatttctgccagagcaagtacagaacgagggccatcgtggcgagaaatcactgtgctccaagctgtattttaccgtaacaagtcattactctgtgcccggagtcatttctcaggatcttaa