



Reference: PTXBC132763

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDC-3; MDC3; disintegrin and metalloproteinase domain-containing protein 23; metalloproteinase-like, disintegrin-like, and cysteine-rich protein 3; ADAM metallopeptidase domain 23		

For research use only.

GENE DESCRIPTION	ADAM metallopeptidase domain 23
GENE ID	8745
ORIGIN SPECIES	Human
ACCESSIONS	BC132763
VECTOR	pDONR223
NCBI SYMBOL	ADAM23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]