

FLJ44006-hypothetical FLJ44006 Gene

Reference: PTXBC132790

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical FLJ44006
GENE ID	400997
ORIGIN SPECIES	Human
ACCESSIONS	BC132790
VECTOR	pDONR223
NCBI SYMBOL	FLJ44006
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaaagcctcaggtgcactaagggtgccaacagaacgatgcacatcagagttgacagtcagactcaaataggacgtgaaggcaacatgtgggtgagagtagcagctatgggtgtaaatatgatcaaatatggggaggtgtccagattcttgacatgggtacaatcagttatttt
tctctcttagtgaacgagtttttgggttttcaatactgctaattttacaggcaattcactacgttccctgggaggtggagtgcccttcctccctgctgtgcgtgggttacacagcctgcctcacttccctgtgggtccttcacacctgcacatgattccattatttgagctcatggcaggaaata
gaacctga