

RFPL4B-ret finger protein-like 4B Gene

Reference: PTXBC132796

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNF211; ret finger protein-like 4B; RING finger protein 211; ret finger protein like 4B		
For research use only.			

GENE DESCRIPTION	ret finger protein-like 4B
GENE ID	442247
ORIGIN SPECIES	Human
ACCESSIONS	BC132796
VECTOR	pDONR223
NCBI SYMBOL	RFPL4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaacgcctgcaagcagagttgtctgtccagtttgcctggatttttctctgttccatttctctctgtacacacgtgttctgtttgattgcacagaggtatatactagaaaaccatgatttagagcgtgtgcccttgtgtcgagacgtgtgaaggtacgtgttgaagaatg
gcaagtgtgagcgtcttaacacttatgaccaagcagcacaatagccgacttgagcaaaagtctgcacgtgagggaggagctccggcattttcgggagagatgtgacctggatgcagccactgccagctccctctgtcttctccaatgatctaagaagcgtcagtgtaagaatccaccac
gatctgacaaaagatcccaggctggcctgtgtctgtggtactccctgtcttctcctccggccaacttactggaggttgaaagtggagaggtgaagtcagtgccctggcgtctgcaaggagccgctgacagaaaagagcaatgatttattccctgagcatggcttctggtacagcatgaa
ggcaggagcaatccatgctaaccacccctggagagaattctgcaagccctgccttgcctgtgtgggaatttctggtgctgacttagaagaatccagtttttgatgttgacaataatgtcctcatatatacatgatgttttctcttggagctttgttccattcttctgtgtgagct
cttgggagaaggaggaggtgcaacgtctgacatctgccatga