

FLJ41649-hypothetical LOC401260 Gene

Reference: PTXBC132803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	long intergenic non-protein coding RNA 951		
For research use only.			

GENE DESCRIPTION	hypothetical LOC401260
GENE ID	401260
ORIGIN SPECIES	Human
ACCESSIONS	BC132803
VECTOR	pDONR223
NCBI SYMBOL	FLJ41649
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcagggaactgggaagctccctccacccatctgatttcagacctgcctcctcagtgcccatgtgtgaccagctatatttagcaagctcaattcagtaactagagcaactggctagcgcctccctccaaggactctgaaggggcagtttctcacaccttttaacatcttccatt
tctccaggcaaccaacacaaaaggcctggaccaggtgggtgagatgctgggctcctacacagaggtctctatttctgttttcagctctccaggccagataagaaaactgagacacagaaaagtcacacaagattatcgagggtcaacacagagaacagtggaaacctgcacctgacccag
gcagcgggaacctcatgaccacaccacaccacaggggaacagcaggggaaaaaatgggtctagcaggacaaaatgacatttccagagtttctgtctgtttggatctgctgcaagtttgcagaagaggcatttccatgagctctggacagtgccagagggcagag
gtgtctccactgcaagtgctggtggtggtccacgggctcctgccaatgtgaggtttgcagcaactgtatataccatggtctgtcaccagcttctgaattcggggcagctggagccacactagtggtccttctgttaggggcagcagctgtctcctggtctcagaatcacagctgtggggtg
caaccttaagccaaaagcctag