

FLJ41649-hypothetical LOC401260 Gene

Reference: PTXBC132805

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC401260
GENE ID	401260
ORIGIN SPECIES	Human
ACCESSIONS	BC132805
VECTOR	pDONR223
NCBI SYMBOL	FLJ41649
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcagggaactgggaagctccctccaccccatctgatttcagacgtgcctcctcagtgcccatgtgtgaccagctatatttagcaagctcaattcagtaactagagcaactggctagcgcctccctccaaggactctgaaggggcagtttctcacaccttttaacatcttcatt
tctccaggcaaccaacacaaaaggcctggaccaggtgggtgagatgctggggtcctacacagaggtctctatttctgttttcagctctccaggccagataagaaaactgagacacagaaaagtcacacaagattatcagggtcaacacagagaacagtggaaacctgcacctgacccag
gcagcgggaacctcatgaccacaccacaccacaggggaacagcaggggaaaaaatgggtctagcaggacaaaatgacattccagagtttctctgtgttggatctgctgcagtttcagttttgccaagaaggcatttccatgagctctggacagtgccagagggcagag
gtgtccactgcagtgctggtggtggtccacgggctcctgccaatgtgaggtttgcagcaactgtatataccatggtctgtcaccagcttctgaattcggggcagctggagccacactagtggtccttctgcaggggcagcagctgtctcctgtctcagaatcacagctgtgggatg
caaccttaagccaaaagcctag