

C20orf191-chromosome 20 open reading frame 191 Gene

Reference: PTXBC132806

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 20 open reading frame 191
GENE ID	149934
ORIGIN SPECIES	Human
ACCESSIONS	BC132806
VECTOR	pDONR223
NCBI SYMBOL	C20orf191
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaagttcaggtaccctcccaacgaaggagcattcagcacagaacaaagtcattatcctcctcactctgtaaagtatacgtttccagcaccaccaccagcaggatccagcatttgaggcaaacatgaagctccatcctccaattctggggcaaccgtgtggagatgatcaaatg
cttcaccttcaaaacttcaagggaagagttatagagtgtagtcgtgtagatcgagaaattgcaaaagtagaacagcagatccttaactgaaaaagaacaagtaaaagtcttcgtctaa