

FLJ41170-hypothetical LOC440200 Gene

Reference: PTXBC132807

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC440200
GENE ID	440200
ORIGIN SPECIES	Human
ACCESSIONS	BC132807
VECTOR	pDONR223
NCBI SYMBOL	FLJ41170
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaaacaaagggttagggatcccgatgccaggtttcagtgctggtgtcctcacaactggctatccctatgccctgctgtcctcagtggtggacactggacctggactgacccctgggacaggaggattcaaggtgtctgtctcttttgattcttttctcttctgccaggaa
agatactgatctctgttcttggtaagtccaagaaccatcctaagttctgtgccccctcagctgtaaaagggaagtacatttcattttttattctgtaaaactctcgggtgtgtgccaatggccatgactgatgatgagcacatgtgtgccccctgccccctggagcgcacatggtctccagc
cagagaccgcgtgggagaaatcaggggttcactcctggtcggagggtgagcatctgcctctgcatgcaggaaggcatctcatgaaccccaaggcctggcagccctgcacatggaaggagtcactctctccatgtgggtgagccacgctggccttggcattcacgtgttccctc
cacctgtctccagcgtga