



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LIPC-lipase, hepatic Gene

Reference: PTXBC132825

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HDLQC12; HTGL; LIPH; hepatic triacylglycerol lipase; Triacylglycerol lipase; hepatic lipase; lipase member C; lipase, hepatic; lipase C, hepatic type		

For research use only.

GENE DESCRIPTION	lipase, hepatic
GENE ID	3990
ORIGIN SPECIES	Human
ACCESSIONS	BC132825
VECTOR	pDONR223
NCBI SYMBOL	LIPC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr