

LOC284100-hypothetical protein LOC284100 Gene

Reference: PTXBC132826

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC284100
GENE ID	284100
ORIGIN SPECIES	Human
ACCESSIONS	BC132826
VECTOR	pDONR223
NCBI SYMBOL	LOC284100
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaccaattgtttacagctctgtccacatgatcattgtaataaggttgagactgagctaaagttaactctggcgatgttctggatgcactggacaaacacctcattccagcagctaccactggaagggtctaccacaggtatctggcagaatttggcaggaatgacaggaagga
ggtaggcagagaacagcttggtggttacaagctgtagtgatattgcaacgggtttctccatgttgccaggctgggctcgaactcctggttcaagtga