



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LOC400707-hypothetical gene supported by AK124070 Gene

Reference: PTXBC132839

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

For research use only.

GENE DESCRIPTION	hypothetical gene supported by AK124070
GENE ID	400707
ORIGIN SPECIES	Human
ACCESSIONS	BC132839
VECTOR	pDONR223
NCBI SYMBOL	LOC400707
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]