

C8orf74-chromosome 8 open reading frame 74 Gene

Reference: PTXBC132852

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C8orf74; chromosome 8 open reading frame 74		
For research use only.			

GENE DESCRIPTION	chromosome 8 open reading frame 74
GENE ID	203076
ORIGIN SPECIES	Human
ACCESSIONS	BC132852
VECTOR	pDONR223
NCBI SYMBOL	C8orf74
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcactctaacacccagggagtgaaagaagtcttcaattcagagaccacaaggtcgggagcgctcgggaggtctgaactgggagagttgacgaacagagagactccggaggagcatcctgtggacaccctctacgagagcatcttgcagtgaggcaaggtt
cccatgggtgagggtgcccaggtggtcaagttcacagaagagctgtaagggaacaaaggctgtccattactgaggtctgtgacgatcctggggaacagcttagagattaccggggccatttcaacaccacccacctgtggccctctgtgactacttccaccacaccttcatccg
cactacaaactaccagtatgtcctgggccaaggaccagcaggtcgacctgacctgtgcccacctggaggtgtgcatgccacccatccctcccgctggccgaggccatggacaggagactgtgtgatccacgagcagcaggtggccacctgacggaggccgagccacagaagcg
cgccgacgtgctcctgaaaaggcgtgctgctggagcgggagaactcgtgcagaaggcgtgctgctgcccgcgcctgctgcagccggccaggtcctgagagacaggaggtggagagcctcatctgccaggcagttccacccagatggagctcctgaggagctgctg
cagcggcagatccagaacacattgcccattggacctgaagcttcagaagaagactctgaacctcaacgccccaccctatccgcccccatcaccagccacgcaggccaggaggaagccctgaagccccaaagagcgagcaaggaaagcgaaggcaaggaaatag