

C16orf86-chromosome 16 open reading frame 86 Gene

Reference: PTXBC132878

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C16orf86; chromosome 16 open reading frame 86		
For research use only.			

GENE DESCRIPTION	chromosome 16 open reading frame 86
GENE ID	388284
ORIGIN SPECIES	Human
ACCESSIONS	BC132878
VECTOR	pDONR223
NCBI SYMBOL	C16orf86
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtccaaagagaaagcctgtcaagtgggggctctcgggggaaggagctgcagcccggggtggggcacacaccacagtaggattcctaaatcctgcactcctgcagggaggattctcgtcccgagccactaccatcactctcctccgttctcaagagccccgttaggga
gcagcaggagagtgaggagcttcagggtcttaggggctgggctcactggccttaagcctctgggctggggtctcattgaccaggaaacgggacaggcgtcagccctgggctgtgggagggctgggctattccctggccagggaatgggagccagggtggagcctggctc
aagtctcttgcctggtcaggtctcagcctccgggcttctgcccatttaaggctgaagctgagctgccacccaagctgccgctgcaggagaggagccagaggacagccagagtgagccctcaccatctgccaacagcacaaaaagccaagaagcgcaagagcctgg
gggtcccgctgctccagctgtggccagcatggtgtctgcaccttagagacattgaggtggagcgtgagtgaggcagtcctggactgttcttctccctgctgtggggcccgacatggcacaacctggcgctctag