

YPEL2-yippee-like 2 (Drosophila) Gene

Reference: PTXBC132884

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FKSG4; protein yippee-like 2; DiGeorge syndrome-related protein; yippee like 2		

For research use only.

GENE DESCRIPTION	yippee-like 2 (Drosophila)
GENE ID	388403
ORIGIN SPECIES	Human
ACCESSIONS	BC132884
VECTOR	pDONR223
NCBI SYMBOL	YPEL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaagatgacaagatcgaagacttccaggcatatctgccctctgccaccggactacagctgcattcaactcagagctcaactggccaatcatgatgaactaattccaagtcattccaaggaagtcaaggacgagcatacctttaaactcagtagttaatgtgggctgtggcctgc
agaagagcagatgttctaacaggactgcagctgcagacatttacttgaaaactgcaaaaccactctggctggaaatacgaacatctttgaaagcagccagaataataaagaaggcaatacatcattgaactagcacacatgatcaaggacaatggctgggactga