



Reference: PTXBC132889

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACH2; CTE-1; LACH2; acyl-coenzyme A thioesterase 1; CTE-I; CTE-Ib; inducible cytosolic acyl-coenzyme A thioester hydrolase; long chain acyl-CoA hydrolase; long chain acyl-CoA thioester hydrolase; acyl-CoA thioesterase 1		

For research use only.

GENE DESCRIPTION	acyl-CoA thioesterase 1
GENE ID	641371
ORIGIN SPECIES	Human
ACCESSIONS	BC132889
VECTOR	pDONR223
NCBI SYMBOL	ACOT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]