



Reference: PTXBC132893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	V-set and immunoglobulin domain-containing protein 8; V-set and immunoglobulin domain containing 8		

For research use only.

GENE DESCRIPTION	V-set and immunoglobulin domain containing 8
GENE ID	391123
ORIGIN SPECIES	Human
ACCESSIONS	BC132893
VECTOR	pDONR223
NCBI SYMBOL	VSIG8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]