

LOC399900-hypothetical gene supported by AK093779 Gene

Reference: PTXBC132894

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical gene supported by AK093779
GENE ID	399900
ORIGIN SPECIES	Human
ACCESSIONS	BC132894
VECTOR	pDONR223
NCBI SYMBOL	LOC399900
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattgggaaggtggaagggtgagacgctccttaggtgtcgcagcgcatcttctcctcagagtcaggaaggagcgcgatggagaagtggggcgacggcggtttcgggtcgcgccctgcccttcgactcctcgaccacgtctgcctttggctgaaagcaaatatggctcc
gtgtgcgtttctgcaacaaggctggcgccgggtttccctagggatgggacgcctcccgggagtcgtgggggtgagccaccgggctccgagagaatcgtgtgtcgtctcgacccaaggacacatctcgggttagaaaggagaaacgacgtggatctaaaggc
gaagcccatgctgaggactttcaaaatggcctcttatcgccccactcacaggagcttcagggttcgggtggccctgtaggacaaccctgggtgtcattggagagaagctgacgggttaattggcggtgacaggtga