

FLJ42842-hypothetical FLJ42842 Gene

Reference: PTXBC132910

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical FLJ42842
GENE ID	440446
ORIGIN SPECIES	Human
ACCESSIONS	BC132910
VECTOR	pDONR223
NCBI SYMBOL	FLJ42842
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctcagtggtggcgtccatcattcctggtgtgtgtgagagtgttccagaggagactcacatgtgagccagcgggctgagaaggagaccgttctcagtggtgagtggtgactgtccaatcagctcaaggcaggctggggacaaacaggcagaagaaggaggattctct
cccaccttctggagcaggatgccttttctccttgacatcagactacagggtcttggcttttgattctaggactgtacaaatggcctccggggccctcagccctcagcctccaacgaaggctgtgtgtgtggtccctgattctgaggctctggactggactgaggcatgctacggg
cttctctgattctccagcttgtggatggcctatcatgggacttctccactctgtatcacaaaggccaatgccccctaataccttctttcatatatactactggttctgtctgcctggggaaccttgactaatacagatatggagcatttgaatga