

MORN5-MORN repeat containing 5 Gene

Reference: PTXBC133690

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf113; C9orf18; MORN repeat-containing protein 5; MORN repeat containing 5		
For research use only.			

GENE DESCRIPTION	MORN repeat containing 5
GENE ID	254956
ORIGIN SPECIES	Human
ACCESSIONS	BC133690
VECTOR	pDONR223
NCBI SYMBOL	MORN5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtacacaggagcaaatatcggggaatatgtagatgggaggatggaggccaagccaagtacatcctcctacgaaacaatatgttggggaatgaaggatggcatgtttcacggcgagggaaccctgtacttcccagcgggaagccaatacagccatttgggaaaa
cggattggccataaaggacacatatacgtttcagatgggctgcactatgatgagaaaaactggcattactgcgacggctatgatcggagggtttacacagagatcctcaatggcttgaagcctgcaggatggctcaactcaccataatggaccacctagaaaaatcccaagggtatta
cgatttggagacggcttctataaccagtcacgagggtatgtaaggactataggaaccgcttcttaagaaacgcagatgatgacgagcatgagtgatcacccgtacctgtcgaagggttag