

TESSP5-testis serine protease 5 Gene

Reference: PTXBC133695

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TESSP5; serine protease 45; testis serine protease 5; protease, serine 45		
For research use only.			

GENE DESCRIPTION	testis serine protease 5
GENE ID	377047
ORIGIN SPECIES	Human
ACCESSIONS	BC133695
VECTOR	pDONR223
NCBI SYMBOL	TESSP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccgccattggccctgggaggtgagcctcggatggaatgagcacgtgtgtgggggcccctcattgacccagctgggtgtgactgcggcccaactgcacccaaggaccaaagagctcagtggtgctggcacctcaagctgcagcccatgaactcagcagggccctctgggtccctgtgaggacattatgcacccaagctggtggccggccctcatcatgggtgacgtgtccctgtccacctcaaacacctgcaccttcagctccacagccaactccatgctgacccagagctgcaggagctgagggtgttatcatggacaacaagaggtgtgaccggcattacaagaagtccttctccccagttgtccctgtctctggggacatgatctgtgcccaattatgggaaaactgtgctatggggtctgagggccattggcttgaagttgagggcagatggattctggctgggggtgtgtctctgggaaaaggcctgcgtgaaggcaca gaatccaggtgtgtacaccgcacccaataacacaaatggatcaagaagcaaatgagcaatggagccttctcaggtccctgtgctcctgcctcctgttctctgtgctgctgcagccccagatgggtcctga