



Reference: PTXBC133707

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACCN2; ASIC; BNaC2; acid-sensing ion channel 1; Cation channel, amiloride-sensitive, neuronal, 2; acid sensing (proton gated) ion channel 1; acid-sensing ion channel 1a protein; amiloride-sensitive cation channel 2, neuronal; brain sodium channel 2; acid		

For research use only.

GENE DESCRIPTION	amiloride-sensitive cation channel 2, neuronal
GENE ID	41
ORIGIN SPECIES	Human
ACCESSIONS	BC133707
VECTOR	pDONR223
NCBI SYMBOL	ACCN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]