

ADAMTSL5-ADAMTS-like 5 Gene

Reference: PTXBC139833

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	THSD6; ADAMTS-like protein 5; ADAMTSL-5; thrombospondin type-1 domain-containing protein 6; thrombospondin, type I, domain containing 6; ADAMTS like 5		

For research use only.

GENE DESCRIPTION	ADAMTS-like 5
GENE ID	339366
ORIGIN SPECIES	Human
ACCESSIONS	BC139833
VECTOR	pDONR223
NCBI SYMBOL	ADAMTSL5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggggcgtatggcgctacgtgcttaattggcactgggtgtcagccaccaggacgtacagggcgccggcagcgtatgtgtctacaccgagacagggcccccaggagacattgcaagcagccggggccacccctccatgacctgctcctacaggtcctcctcaggagc
ccaaccctggcagcagtttgggtctgctccctcgggagcgtacagccctccaggctcgtgtcagggccctgggctggccctgaggcagcctcagccccgggggtggagcctcagcccccgagccctgctgtacccctgcacagaccccaacgtggccccagacc
cctgccaccctgcccagaccccgccggccgccccaccgactactcactattgcggcagtgactttgtgtccagggcccgagtgctggggccaccaccaccagggccagggagaccgctatgaggtgctgcatccagctcgtctacaagaaccctcgcctcctcgggcacgcgagt
acgtgtggggccagggccactgcccctgcccgtgctgacccaccgggactacctgatggctgtccagcgtctgtcagccccgacggcacacaggaccagctgtgtgtccccacggcctacgcccggccctggagccctcggaggacagccgcatacgcctgactgcc
cggcgtgtctcgtgctga