

FLJ45422-FLJ45422 protein Gene

Reference: PTXBC139904

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FLJ45422 protein
GENE ID	441140
ORIGIN SPECIES	Human
ACCESSIONS	BC139904
VECTOR	pDONR223
NCBI SYMBOL	FLJ45422
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtggctccgacagacacctggagaacgggaaggagcgtgacagcagcggatccccaaaggcacatgtgacccagcacccatctctgaccatgaggccacccctgaggtgctgggcccctggcctctaccctgaggagatcacactgacctggcagcaggatggggagg
accagacccaggacacggagctgtggagaccagcctgcaggggacggaacctccagaagtgggtgctgtagtggcttcggagaggagcagagatacatgtgccatgtgcagcatgaggggctgccagagccccctaccctgagatgggagccgtctctcagccacc
atccccatcgtggcgtctgtgctggcctgtttctcctggagctgtggtcactggagctgtggtgctgctgcgatgtggaggaagaaagctcagcagcaattgtgctcagctctgatgcctcatgatactgtaaaggagactatgcctgttctgttctgtctgctgtctct
ttctccctgggtgtccctcatctctgacagcagcaggagtcattttctgtcatttaacccacaagggtgaagcagccccctgcacacagaagtctgtgtattaa