



## Reference: PTXBC139907

|                               |                                          |                     |                |
|-------------------------------|------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                      | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | DIPOA; EDM5; HOA; OADIP; OS2; matrilin-3 |                     |                |
| <b>For research use only.</b> |                                          |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | matrilin 3                                                                                       |
| <b>GENE ID</b>          | 4148                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC139907                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | MATN3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]