

SLIT1-slit homolog 1 (Drosophila) Gene

Reference: PTXBC139909

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MEGF4; SLIL1; SLIT-1; SLIT3; slit homolog 1 protein; multiple EGF-like domains protein 4; multiple epidermal growth factor-like domains protein 4; slit homolog 1; slit guidance ligand 1		

For research use only.

GENE DESCRIPTION	slit homolog 1 (Drosophila)
GENE ID	6585
ORIGIN SPECIES	Human
ACCESSIONS	BC139909
VECTOR	pDONR223
NCBI SYMBOL	SLIT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctgactccgggtgggggtcctcggggggccggtccggccggagctctgctgctgtgggcagccgctggcgctgggtgctcggcggtgccccgccctctgcacctgcaccggaaccacgggtgactgccacggcaggggctgcaggccattccaagaata
 tactctggaacaccgagcgctgtga