

RHD-Rh blood group, D antigen Gene

Reference: PTXBC139922

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	blood group protein RHD; RhD polypeptide; RhD blood group antigen; CD240D; IIIc; RH30; RHCED; RHDVA(TT); RHDel; RHPH; RHXIII; Rh4; RhDCw; RhII; RhK562-II; RhPI; blood group Rh(D) polypeptide; D antigen (DCS); RH polypeptide 2; Rh blood group CcEe antige		

For research use only.

GENE DESCRIPTION	Rh blood group, D antigen
GENE ID	6007
ORIGIN SPECIES	Human
ACCESSIONS	BC139922
VECTOR	pDONR223
NCBI SYMBOL	RHD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgagctctaagtaccgcggtctgtccggcgtcctgcccctctgggcccctaaactggaagcagctctcattctctcttcttattttaccctatgacgcttcttagaggatcaaaaggggctgtggcatcctcaagtggccaagatctgaccgtgatggcggccattggcttgg
gcttctcaccctcagttfccggagacacagctggagcagtggtgcccctcaacctctcatgctggcgttgggtgagtgaggcaatcctgctggagcgtctctgagccagttccctctgggaaggtggtcatcacactgttcagtattcggctggccaccatgagtccttctgctgctg
atcagtgatgctgtcttgggaaggtcaactggcgagttgggtgatggtgctggtggagtgacagctttaggcaacctgaggatggtcagtaatacttcaacacagactaccacatgaacatgatcacatctacgtttcgcagcctatttggctgtctggtgctggtgc
ctgccaagcctctaccggagggaacggagataaagatcagacagcaacgataccagttgtctccatgctggcgccctctctgtggatgtctgccaagttcaactctgctgctgagaagtccaatgaaaggaagaatgccgtgttcaacactactatgctgtagcagta
gctgtgtgacagccatcaggggtcatcctggctcaccocccaagggaagatcagcaagactatgtgacagtgcggtgtggcaggagcgtgctgtgggtacctgtgtaacctgacccctctccgtggcttccatggtgctgggtctgtggtggctgatctcctcgggga
gccaagtacctgcccgttctcattggctgttggatftaa
```