

## Human UHRF1(400-660) Recombinant Protein

PX-P1062-10

### OVERVIEW

<b>SIZE</b>	10 ug
<b>ORIGIN SPECIES</b>	Human
<b>FRAGMENT</b>	Partial*
<b>PROTEIN DELIVERED WITH TAG</b>	Yes
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	55,97 kDa
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	E.coli
<b>PURITY</b>	0,9
<b>PROTEIN ACCESSION</b>	AAF28469.1
<b>FORM</b>	liquid
<b>BUFFER</b>	TrisHC 50mMl, 10mM reduced glutathion pH8 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>STABILITY &amp; STORAGE</b>	

### PROTEIN SEQUENCE

MSPILGYWKIKGLVQPTRLLEYLEEKYEEHL YERDEGDKWRNKKFELGLEFPNLPYYIDGDVKTQSMAIRYIADKHN  
MLGGCPKERAESMLEGAVLDIRYGVSR IAYSKDFETLKVDFLSKLP EMLKMFEDRLCHKTYLNGDHVTHPDFMLYDALD  
VVLYMDPMCLDAFPKLVCFKKRIEAI PQIDKYLKSSKYIAWPLQGWQATFGGGDHPPKSDLEVL FQGPLGSPEFMACVGR  
TKECTIVPSNHYGPIPGIPVGT MWRFVQVSESGVHRPHVAGIHGRSNDGAYSLV LAGGYEDD VDHGNFFTYTGSGGRDL  
SGNKR TAEQSCDQKLTNTNRALALNCFAPINDQEGAEAKDWRSGKPVRRVVRNVKGGKNSKYAPAEGNRYDGIYKVVKYWP  
EKGKSGFLVWRYLLRRDDDEPGPWTKEGKDRIKLG LTMQYPEGYLEALANREREKENS KREEEEQQEGGFASPRTGK GK  
WKRKSAGGGPSRAG

### More info

#### SPECIFICATION

#### GENE ID

29128

#### SWISSPROTID

Q96T88

#### NCBI GENE ALIASES

hUHRF1, ICBP90,  
RNF106, Np95, huNp95,  
hNP95

#### UNIPROT ID

Q96T88

#### UNIPROT LINK

<http://www.uniprot.org/uni/prot/Q96T88>

#### SYNONYMS

ICBP[400-660],  
transcription factor ICBP90