

Toxoplasma gondii GRA2(21-185)(76K) Recombinant Protein

PX-P1094-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Toxoplasma gondii
FRAGMENT	Partial
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	25,78 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,8
PROTEIN ACCESSION	XP_002366395.1
FORM	liquid
BUFFER	PBS, imidazole 300mM 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MHHHHHSSGLVPRGSGMKETAALKFERQHMDSPDLGTDDDDKSPGFSSTMGMRAAEFSGVVNQGPVDVPFSGKPLDERA
VGGKGEHTPPLPDERQQEPEEPVSQRASRVAEQLFRKFLKFAENVGQHSEKAFKKAKVVAEKGFTAAKTHTVRGFKVAKE
AAGRMVTVGKKLANVESDRSTTTTQAPDSPNGLAETQAPVEPQQRAAHVVPDFSQSDPNSSSVDKLAAALEHHHHHHH

More info

SPECIFICATION

GENE ID

7897424

SWISSPROTID

S8EXV7

NCBI GENE ALIASES

UNIPROT ID

B6KEX2

UNIPROT LINK

<http://www.uniprot.org/uni-prot/B6KEX2>

SYNONYMS

GRA2(21-185)(76K), 28kd antigen