

Toxoplasma gondii GRA2(21-185)(RH) Recombinant Protein

PX-P1095-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Toxoplasma gondii
FRAGMENT	Partial
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	25,86 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,9
PROTEIN ACCESSION	AAA30138.1
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MHHHHHSSGLVPRGSGMKETA AAKFERQHMDSPDLGTDDDDKSPGFSSTMDMRAAEFSGVVNQGPVDVPFSGKPLDERA
VGGKGEHTPPLPDERQQEPEEPVSQRASRVAEQLFRKFLKFAENVGHHSEKAFKKAKVVAEKGFTA AKTHTVTRGFKVAKE
AAGRMVTVGKKLANVESDRSTTTTQAPDSPNGLAETEVPVEPQRAAHVPVPDFSQADPNSSSVDKLAAALEHHHHHHH

More info

SPECIFICATION

GENE ID

SWISSPROTID

P13404

NCBI GENE ALIASES

UNIPROT ID

P13404

UNIPROT LINK

<http://www.uniprot.org/uni/prot/P13404>

SYNONYMS

GRA2(21-185)(RH), 28kd antigen