

## Toxoplasma gondii GRA2(Nt-Ct) (RH) Recombinant Protein

PX-P1096-10

### OVERVIEW

<b>SIZE</b>	10 ug
<b>ORIGIN SPECIES</b>	Toxoplasma gondii
<b>FRAGMENT</b>	Partial
<b>PROTEIN DELIVERED WITH TAG</b>	Yes
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	18,41 kDa
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	E.coli
<b>PURITY</b>	0,9
<b>PROTEIN ACCESSION</b>	XP_002366395.1
<b>FORM</b>	liquid
<b>BUFFER</b>	PBS, imidazole 300mM
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

### PROTEIN SEQUENCE

MHHHHHSSGLVPRGSGMKETAALKFERQHMDSPDLGTDDDDKSPGFSSTMGRAAEFSGVVNQGPVDVPFSGKPLDERAV  
GGKGEHTPPLPDERQPEEASMTVTGKLANVESDRSTTTTQAPDSPNGLAETVPVPEPQRAAHVPVPDFSQNSSVVDK  
LAAALEHHHHHHH

### More info

#### SPECIFICATION

#### GENE ID

7897424

#### SWISSPROTID

S8EXV7

#### NCBI GENE ALIASES

#### UNIPROT ID

B6KEX2

#### UNIPROT LINK

<http://www.uniprot.org/uni-prot/B6KEX2>

#### SYNONYMS

GRA2(Nt-Ct) (RH), 28kd antigen