

Toxoplasma gondii GRA6Ct(174-230)(RH) Recombinant Protein

PX-P1100-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Toxoplasma gondii
FRAGMENT	Partial
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	11,54 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,9
PROTEIN ACCESSION	BAE86889.1
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MHHHHHSSGLVPRGSGMKETAALKFERQHMDSPDLGTTGRRSPQEPSGDGGGNDAGNNAGNGGNEGRGYGGRGEGGAED
DRRPLHPERVNVFDYKLAALAEHHHHHHH

More info

SPECIFICATION

GENE ID**SWISSPROTID**

Q25C69

NCBI GENE ALIASES**UNIPROT ID**

Q25C69

UNIPROT LINK

<http://www.uniprot.org/uni/prot/Q25C69>

SYNONYMS

GRA6Ct(174-230)(RH),
dense granule antigen