

Trichinella spiralis GST24 Recombinant Protein

PX-P1104-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Trichinella spiralis
FRAGMENT	Full
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	25,32 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	>95%
PROTEIN ACCESSION	ABA42914.1
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MAHNHRHKHKLMPLYKLVYFPIRGLAEPIRLLLHDQRVEFLDNRIQKQDWPEIKSQMLFGQVPCLYEDDQPIVQSGAIMR
HLGRRFGLYGNAEEMTYVDQIYEGVVDLRLKYARLIYSDSFHESKGGKFINEVLPDELAKFEKILTGKKYILDDEITFADY
ALAELLDVLLILSSSCLNFALTIIYHSRFMNRPNLKRYLSSDIRKNAKINGNENK

More info

SPECIFICATION

GENE ID

SWISSPROTID

Q0H8U8

NCBI GENE ALIASES

UNIPROT ID

Q0H8U8

UNIPROT LINK

<http://www.uniprot.org/uni/prot/Q0H8U8>

SYNONYMS

GST24, glutathione-S-transferase