

Drosophila SO Recombinant Protein

PX-P1154-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Drosophila
FRAGMENT	Partial (Nter)*
PROTEIN DELIVERED WITH TAG	No
MOLECULAR WEIGHT WITH TAG IF ANY	8,60 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,6
PROTEIN ACCESSION	NP_476733.1
FORM	liquid
BUFFER	Tris-HCl 50mM, NaCl 150mM, EDTA 1mM, DTT 1mM, pH7
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

GPLGSM LQHPATDFYDLAAANAAAVLTARHTPPYSPTGLSGSVALHNNNNNSSTSNNNNSTLDIMAHNGGGAGGGLHLN
SSSNGGGG

More info

SPECIFICATION

GENE ID

35662

SWISSPROTID

Q27350

NCBI GENE ALIASES

ami, med, old, mda, Mdu,
So, Drl, SO, somda,
Dmel\CG11121, CG11121

UNIPROT ID

Q27350

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q27350>

SYNONYMS

SO, sine oculis