

Cow SAA3 Recombinant Protein

PX-P2056-10

OVERVIEW

SIZE	10 uG
ORIGIN SPECIES	Cow
FRAGMENT	Partial (no SIP)*
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	14,19 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,9
PROTEIN ACCESSION	NP_851359.2
FORM	liquid
BUFFER	PBS, imidazole 400mM, pH7.4
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MRGSHHHHHHGSQRWGTFLKEAGQGAKDMWRAYQDMKEANYRGADKYFHARGNYDAARRGPGGAWAAKVISNARETIQGI
TDPLFKGMTRDQVREDSKADQFANEWGRSGKDPNHFRPAGLPDKY

More info

SPECIFICATION

GENE ID

281474

SWISSPROTID

Q8SQ28

NCBI GENE ALIASES

UNIPROT ID

Q8SQ28

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q8SQ28>

SYNONYMS

SAA3, serum amyloid A 3 precursor