

Human IGFBP1 Nter Recombinant Protein

PX-P2078-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	17.59 kD
DELIVERY CONDITION	The product is shipped at ambient temperature

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	85%
PROTEIN ACCESSION	ACO37638.1
FORM	Lyophilized
BUFFER	PBS, pH 7.5 with 0.2mM β -Met 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MKYLELLPTAAAGLLLLAAQPAMA[^]MDAPWQCAPCSAEKLALCPPVSASCSEVTRSAGCGCCPMCALPLGAACGVATARCARGLSICALPGEQQPLHALTRGQGA
CVQESDASAPHAEEAGSPESPESTEITEEELLDNFHLMAPSEEDHSILWDAISTYDGSKLEHHHHHHH

More info

SPECIFICATION

GENE ID

3484

SWISSPROTID

P08833

NCBI GENE ALIASES

hIGFBP-1, AFBP, IBP1,
PP12, IGF-BP25

UNIPROT ID

P08833

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P08833>

SYNONYMS

IGFBP1 Nter, insulin-like
growth factor binding
protein 1, Insulin-like
growth factor-binding