

Human RNase I Mutant type Recombinant Protein

PX-P2089-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	15.75 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	90%
PROTEIN ACCESSION	NP_002924
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MGKESRAKKFQRQVMDSDSPSSSTYCNQMRRRNMTQGRCKPVNTFVHEPLVDVQNVCFQEKVTCKNGQNCYKSNSSMHITDCRLTNGSRYPNCAYRTS
PKERHIIVACEGSPYVPVHFDASVEDSTLEHHHHHHH

More info

SPECIFICATION

GENE ID

6035

SWISSPROTID

P07998

NCBI GENE ALIASES

RIB1, RNS1, RAC1

UNIPROT ID

P07998

UNIPROT LINK

<http://www.uniprot.org/uni/prot/P07998>

SYNONYMS

RNase I Mutant type,
ribonuclease, RNase A
family, 1 (pancreatic)