

Streptococcus 6His-Protein G

PX-P2102-10mg

OVERVIEW

SIZE	10 mg
ORIGIN SPECIES	Streptococcus
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	20.96 KDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	90%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MGKGETTEAVDAATAAEKVKFYANDNGVDGEWYDDATKFTVTEKPEVIDASELTPAVTTYKLVINGKTLKGETTEAVDAATAAEKVKFYANDNGVDGEWYDDATKFTVTEKPEVIDASELTPAVTTYKLVINGKTLKGETTTKAVDAETAEKAFKQYANDNGVDGVWYDDATKFTVTEHHHHHH

More info

SPECIFICATION

GENE ID

SWISSPROTID

P19909

NCBI GENE ALIASES

UNIPROT ID

UNIPROT LINK

SYNONYMS

6His-Protein G