

## Human Resistin-like alpha (HIMF) Recombinant Protein

PX-P3011-10

### OVERVIEW

<b>SIZE</b>	10 ug
<b>ORIGIN SPECIES</b>	Human
<b>FRAGMENT</b>	
<b>PROTEIN DELIVERED WITH TAG</b>	Yes
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	36 kDa
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	mammalian
<b>PURITY</b>	85%
<b>PROTEIN ACCESSION</b>	
<b>FORM</b>	Frozen
<b>BUFFER</b>	Tris-HCl 0.1M, Glycine 0.1M, pH8.5 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>STABILITY &amp; STORAGE</b>	

### PROTEIN SEQUENCE

MKTATCSLLICVFLQLMVPVNTDGTLDIIGKKKVKELLAHQDNYPSAVRKTLSCTNVKSMSKWASCPAGMTATGCSGFCAGSWEIQNENICNCLCLIVDWAY  
ARCCQLSLEVLFGQDPDKTHTCPPCPAPPELLGGPSVFLFPPKPKDQLMISRTPVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQNSTYRVVSVLTVLH  
QDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKS  
RWQQGNVFSVLSLHEALHNHYTQKSLSLSPGK

### More info

#### SPECIFICATION

**GENE ID**

**SWISSPROTID**

**NCBI GENE ALIASES**

**UNIPROT ID**

Q9BQ08

**UNIPROT LINK**

<http://www.uniprot.org/uni-prot/Q9BQ08>

#### SYNONYMS

Resistin-like beta; Colon and small intestine-specific cysteine-rich protein; Colon carcinoma-related gene protein; Cysteine-rich secreted protein A12-alpha-like 1; Cysteine-rich secreted protein FIZZ2;