

Chinese sturgeon Luteinizing hormone(Luteinizing) Recombinant Protein

PX-P3017-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Chinese sturgeon
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	40.94 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	90%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MPASMPLLLLLFSALILSHSSSLRLCEPVNETISAEKEECPKCLLIQTSICSGSCTKDPVFKSALSTVQQHVCTYKDLRFVTVLTPDCPPGVDPHFTFPLALSCECSLC
RMESDCTIQSVGSDCMGELAIQNYDDDDKDKTHTCPPAPELLGGPSVFLFPPKPKDQLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPRE
EQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVL
DSDGSEFLYSLKTLVDSKSRWQQGNVFSCSVLHEALHNHYTQKSLSLSPGK

More info

SPECIFICATION

GENE ID

SWISSPROTID

NCBI GENE ALIASES

UNIPROT ID

B2CKR4

UNIPROT LINK

<http://www.uniprot.org/uniprot/B2CKR4>

SYNONYMS

ICSH; Lutropin; Interstitial Cell Stimulating Hormone