

Human BST2 Recombinant Protein

PX-P3037-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	39.49 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	80%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

SEACRDGLRAVMECRNVTHLLQQLTEAQKGFQDVEAQAATCNHTVMALMASLDAEKAQGQKKVEELEGITTLNHLKQDASAEVERLRRENQVLSVRIADK
KYYPSQDSS

More info

SPECIFICATION

GENE ID

SWISSPROTID

NCBI GENE ALIASES

UNIPROT ID

Q10589

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q10589>

SYNONYMS

CD317; Tetherin; HM1.24
antigen; Bone Marrow
Stromal Cell Antigen 2;
BST2