

Human CAV1 Recombinant Protein

PX-P3038-100

OVERVIEW

SIZE	100 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	12.71 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	95%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

SGGKYVDSEGHLYTVPIREQGNIYKPNNKAMADELSEKQVYDAHTKEIDLVNRPKHLNDDVVKIDFEDVIAEPEGTHSFDGIWKASFTTFTVTKYWFYRLLS

More info

SPECIFICATION

GENE ID

SWISSPROTID

NCBI GENE ALIASES

UNIPROT ID

Q03135

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q03135>

SYNONYMS

VIP21; Caveolae Protein
22kDa; Caveolin 1; CAV1