

Human FLOT1 Recombinant Protein

PX-P3042-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	31.78 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	85%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

TLKDIHDDQDYHLHSLGKARTAQVQKDARIGEAEAKRDAGIREAKAKQEKVSAQYLSEIEMAKAQRDYELKKAAYDIEVNTRRAQADLAYQLQVAKTKQQIEEQ
RVQVQVVERAQQVAVQEIEARREKELEARVRKPAEAERYKLERLAEAEKSQLIMQAEAEAAASVRMRGEAEFAIGARARAEAEQMAKKAFAFQLYQEAAQL
DMLLEKLPQVAEEISGPLTSANKITLVSSGSGTGMGAAKVTGEVLDILTRLPEVERLTGVSISQVNHKPLRTA

More info

SPECIFICATION

GENE ID

SWISSPROTID

NCBI GENE ALIASES

UNIPROT ID

O75955

UNIPROT LINK

<http://www.uniprot.org/uni/prot/O75955>

SYNONYMS

Flotillin-1; FLOT1