

Human TSG101 Recombinant Protein

PX-P3051-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	17.44 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	80%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MAVSESQLKMKMVSKEYKRDLTVRETVNVITLYKDLKPVLDYSYVFDNGSSRELMNLTGTIPVPPYRGNTYNIPICLWLLDTYPYNPPICFVKPTSSMTIKTGKHVDAN
GKIYLPYLHEWKHPQSDLLGLIQVMIVVFGDEPPVFSRP

More info

SPECIFICATION

GENE ID

SWISSPROTID

NCBI GENE ALIASES

UNIPROT ID

Q99816

UNIPROT LINK

<http://www.uniprot.org/uni/prot/Q99816>

SYNONYMS

TSG10; VPS23; Tumor

Susceptibility Gene 10;

ESCRT-I complex subunit

TSG101;Tumor

Susceptibility Gene 101