

Listeria monocytogenes Invasion associated protein p60 Recombinant Protein

PX-P3056-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Listeria monocytogenes
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	48.7 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	85%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH 7.5 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MSTVVVEAGDTLWGIAQSKGTTVDAIKKANLLTDDKIVPGQKLQVNEVAAAETKESVSATWLNVRSGAGVDNSIITSIKGGTKVTVESTESNGWNKITYNDGETGFVNGKYLTDKIVASTPVAPTQEVKKETTTQQAAPAAETKTEVKQTTQATTPAPKVAETKETPVVDQNAATHAVKSGDTIWALSVKYGVSVDIMSWNNLSSSIYVGGKLAIKQTANTATPKAEVKKEAPAAEKQAAPVVKENTNTNTNATTEKKETTQQQTAPKAPTEAAKPAPAPSTNTNANKTNTNTNTNTSTPSKNTNTNTNSNTNTNSNTNANQGSSNNNSNSSASAIIEAQKHLGKAYSWGNGPPTTFDCSGYTKYVFAKAGISLPRTSQAQYASTTRISESQAKPGDLVFFDYGGSGISHVGIYVNGQMINAQDNGVKYDNIHSGWGKYLGVGFRVLEHHHHHH

More info

SPECIFICATION

GENE ID

SWISSPROTID

NCBI GENE ALIASES

UNIPROT ID

Q84DP3

UNIPROT LINK

<http://www.uniprot.org/uni/prot/Q84DP3>

SYNONYMS

iap; Invasion associated protein p60