

Rainbow trout rtCD64s Recombinant Protein

PX-P3064-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Rainbow trout
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	34.74 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	90%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS pH 7.5 +urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

QVLPQSAPQTAVAELVKGLPWIFSGESVHLKCSVPGNVSAEWRYRWFRRGGEQLQSESYFVLWKARPPQSGKYYCRGLRRTTWINTQHTLHSLPIEIEVDGGWAIL
QAPPLPMLVGETMTLTCRVRNPNKLTEVILYKDGVELQIRGPELRVTNLTQHGHGSYWCASWDGRRETNSVISVAAPVSIVEVLTEPMLIVPNPLIHKDRML
LVCHVQLNARQPVPHIHYFYQDGLSLGPASSQDKITVLRDSGQYWCKASVPTLGLKRLSKPLGYGRVTGEPSEHLRMRNQKIIVVPLGSHHHHHH

More info

SPECIFICATION

GENE ID
SWISSPROTID
NCBI GENE ALIASES
UNIPROT ID
UNIPROT LINK

SYNONYMS

rtCD64s