

Human CD19 recombinant protein

PX-P4004-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	34.12kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,800000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	150 mM NaCl, 50 mM Tris-HCl, pH 7.5 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MPPPRLLFFLLFLTPMEVRPEEPLVVKVEEGDNAVLQCLKGTSDGPTQQLTWSRESPLKPFLKL
SLGLPGLGIHMRPLAIWLFIFNVSQQMGGFYLCQPGPSEKAWQPGWTVNVEGSGELFRWNV
SDLGGLGCGLKNRSSEGPSKGLMSPKLYVWAKDRPEIWEGEPPCLPPRDSLNSLSQDLT
MAPGSTLWLSGVPDVS SRGPLSWTHVHPKGPKSLLSLELKDDRPARDMWV METGLLLPRA
TAQDAGKYCHRGNTMSFHLEITARPVLWHWLLRTGGWKENLYFQGSHHHHHH

More info

UNIPROT ID

P15391

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P15391>

SYNONYMS

B-lymphocyte antigen
CD19, CD19