

Mus musculus Leukocyte surface antigen CD47 recombinant protein

PX-P4009-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Mus musculus
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	18.85kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,800000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50mM Tris-HCl pH 7.0, 150mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MWPLAAALLGSCCCGSAQLLFSNVNSIEFTSCNETVVIPCIVRNVEAQSTEEMFVKWKLNKS
YIFIYDGNKNSTTTDQNFTSAKISVSDLINGIASLKMMDKRDAMVGNYTCEVTELSREGKTVIEL
KNRTAFNTDQGSACSYYYYEKGCKLVSFSPNEKIGSHHHHHH

More info

UNIPROT ID

Q61735

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q61735>

SYNONYMS

CD47, MER6, OA3,
Antigen Identified By
Monoclonal Antibody 1D8,
Antigenic Surface
Determinant Protein OA3,
Rh-Related Antigen