

Human Alpha-fetoprotein/Beta-actin recombinant protein

PX-P4010-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,9499999999999996
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

More info

UNIPROT ID

P02771

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P02771>

SYNONYMS

A-FP, FETA, HPAPF,
Alpha-Fetoglobulin, Alpha-1-fetoprotein