

Human CD272/B- and T-lymphocyte attenuator(BTLA) recombinant protein

PX-P4013-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	28.73kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,800000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl, pH 8.0, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MKTLPAMLGTGKLFWVFFLIPYLDIWNHIGKESCDVQLYIKRQSEHSILAGDPFELECPVKYCA
NRPHVTWCKLNGTTCVKLEDRQTSWKEEKNISFFILHFEPVLPNDNGSYRCSANFQSNLIESHS
TTLVYTDVKSASERPSKDEMASGSHHHHHH

More info

UNIPROT ID

Q7Z6A9

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q7Z6A9>

SYNONYMS

BTLA, CD272, BTLA1, B- and T-lymphocyte-associated protein,