

Human CD40/ Tumor necrosis factor receptor superfamily member 5(TNFRSF5) recombinant protein

PX-P4015-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	22.32kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,900000000000000002
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl, pH 8.0, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MVRLPLQCVLWGCLLTAVHPEPPTACREKQYLINSQCCSLCQPGQKLVS DCTEFTETEC LPCG
ESEFLDTWNRETHCHQHKYCDPNLGLRVQQGTSETDTICTCEEGWHCTSEACESCVLHRSCS
PGFGVKQIATGVSDTICEPCPVGFFSNVSSAFEKCHPWTSCETKDLVVQQAGTNKTDVVC GPQ
DRLRGSHHHHHH

More info

UNIPROT ID

P25942

UNIPROT LINK

<http://www.uniprot.org/uniprot/P25942>

SYNONYMS

CD40, CDW40, Bp50, P50,
B-cell surface antigen
CD40, CD40L
receptor, TNFRSF5